

DECLARATION OF BLOOD PRESSURE MEASURING DEVICE EQUIVALENCE 2006

A SIGNED COPY WILL BE POSTED ON THE www.dableducational.org WEBSITE

SECTION A - Please complete all items online.

I Tomohiro Kukita Director of Omron Healthcare Europe B.V.
Name of a Company Director Company name

hereby state that there are no differences that will affect blood pressure measuring accuracy between the

Omron BP785 (HEM-7222-Z)
Blood pressure measuring device for which validation is claimed

blood pressure measuring device and the

Omron M6 Comfort (HEM-7000-E)
Existing validated blood pressure measuring device

blood pressure measuring device, which has previously passed the International protocol, the results of which were published as follows

Belghazi J, El Feghali RN, Moussalem T, Rejdych M, Asmar RG
Authors(s)

Validation of four automatic devices for self-measurement of blood pressure according to the International Protocol of the European Society of Hypertension
Title

Vascular Health and Risk Management 2007;3(4):389-400
Publication Year Volume Pages

The only differences between the devices involve the following components:

(When a component is not relevant, both Yes and No should be left blank. Please provide details on any differences below.)

| | | | | |
|---------|----|---|---|--|
| Part I | 1 | Algorithm for Oscillometric Measurements | Yes <input type="checkbox"/> | No <input checked="" type="checkbox"/> |
| | 2 | Algorithm for Auscultatory Measurements | Yes <input type="checkbox"/> | No <input type="checkbox"/> |
| | 3 | Artefact/Error Detection | Yes <input type="checkbox"/> | No <input checked="" type="checkbox"/> |
| | 4 | Microphone(s) | Yes <input type="checkbox"/> | No <input type="checkbox"/> |
| | 5 | Pressure Transducer | Yes <input checked="" type="checkbox"/> | No <input type="checkbox"/> |
| | 6 | Cuff or Bladder | Yes <input type="checkbox"/> | No <input checked="" type="checkbox"/> |
| | 7 | Inflation Mechanism | Yes <input type="checkbox"/> | No <input checked="" type="checkbox"/> |
| | 8 | Deflation Mechanism | Yes <input type="checkbox"/> | No <input checked="" type="checkbox"/> |
| Part II | 9 | Model Name or Number | Yes <input checked="" type="checkbox"/> | No <input type="checkbox"/> |
| | 10 | Casing | Yes <input checked="" type="checkbox"/> | No <input type="checkbox"/> |
| | 11 | Display | Yes <input checked="" type="checkbox"/> | No <input type="checkbox"/> |
| | 12 | Carrying/Mounting Facilities | Yes <input type="checkbox"/> | No <input type="checkbox"/> |
| | 13 | Software other than Algorithm | Yes <input checked="" type="checkbox"/> | No <input type="checkbox"/> |
| | 14 | Memory Capacity/Number of stored measurements | Yes <input checked="" type="checkbox"/> | No <input type="checkbox"/> |
| | 15 | Printing Facilities | Yes <input type="checkbox"/> | No <input type="checkbox"/> |
| | 16 | Communication Facilities | Yes <input type="checkbox"/> | No <input type="checkbox"/> |
| | 17 | Power Supply | Yes <input type="checkbox"/> | No <input checked="" type="checkbox"/> |
| | 18 | Other Facilities | Yes <input type="checkbox"/> | No <input checked="" type="checkbox"/> |

Brief explanation of differences and further relevant details:

5) The pressure sensor is replaced to a piezo electric sensor (NPS) from a capacitive sensor (CPSU), but the accuracy of blood pressure measurement is equivalent between NPS and CPSU.

10) The “user ID selection” switch, the “Morning/Evening Average” switch and the “Memory” switch are added.

11) The symbol for cuff wrapping guide, the indicator for blood pressure level, the morning hypertension symbol, the morning average symbol, the evening average symbol, the week display, the user ID symbol and the TruRead™ symbol are added.





13) The function to guide cuff wrapping, the function to detect morning hypertension, the function to calculate a weekly averages for measurements taken in the morning and evening and the function to take three consecutive measurements (TruRead™) are included.

14) Stores 100 readings each for two users instead of 90 readings for one user.

SECTION B - Complete all items, bar signatures and seal, online and print. Sign and seal it then send the original along with manuals for both devices to our address below.

| | | |
|-----------------------|--|--|
| Signature of Director | <u>Tomohiro Kukita.</u> | Company Stamp/Seal |
| Name | <u>Tomohiro Kukita</u> | |
| Date | <u>8th June, 2011</u> | OMRON HEALTHCARE EUROPE B.V. ' Kruisweg 577 NL-2132 NA Hoofddorp → P.O. Box 2150 NL- 2130 GL Hoofddorp Tel. +31 - 20 354 82 00 Fax +31 - 20 354 82 01 |
| Signature of Witness | <u>J. Meijer-Dul</u> | |
| Name | <u>J. MEIJER-DUL</u> | |
| Address | <u>Omron Healthcare Europe B.V., Kruisweg 577, 2132NA Hoofddorp, The Netherlands</u> | |

Comparison of the Omron BP785 (HEM-7222-Z) with the Omron M6 Comfort (HEM-7000-E)

| Devices | Omron BP785 (HEM-7222-Z) | Omron M6 Comfort (HEM-7000-E) |
|-------------------|---|---|
| Pictures |  |  |
| Display |  |  |
| Validation | | ESH-IP 2002 |
| Device 1 Criteria | <p>Measurement</p> <p><i>Sensors</i></p> <p>Pressure sensor: 2nd sensor for dual check 5</p> <p>Buttons/Switches</p> <p><i>Measurement Records</i></p> <p>User ID 10</p> <p><i>Analysis</i></p> <p>Morning/Evening Average 10</p> <p>Display/Symbols/Indicators</p> <p><i>Preparation</i></p> <p>Correct cuff wrapping indicator 11, 13, 18</p> <p>Multiple measurements (3) 11, 13</p> <p><i>Measurement Records</i></p> <p>User (A or B) + Guest 11</p> | |

| Devices | Omron BP785 (HEM-7222-Z) | Omron M6 Comfort (HEM-7000-E) |
|--------------------------------------|---|---|
| Device 1 Criteria (continued) | <p>Display/Symbols/Indicators (continued)</p> <p><i>Post Measurement</i></p> <p>Morning hypertension (rising sun) 11, 13</p> <p>Daytime weekly average (sun) 11, 13</p> <p>Night-time weekly average (moon) 11, 13</p> <p>Week indicator 11, 13</p> <p><i>Settings</i></p> <p>Sensor cross check (LED) 5, 18</p> <p>Algorithms</p> <p><i>Averages and Differences</i></p> <p>Daytime weekly average × 8 weeks 13</p> <p>Night weekly average × 8 weeks 13</p> <p><i>Parameter Settings</i></p> <p>Correct cuff wrapping detection 13</p> <p>Sensor cross check 5, 18</p> | |
| Same Criteria | <p>Measurement</p> <p><i>Accuracy</i></p> <p>BP accuracy ± 3 mmHg or ± 2%^{Query 1} 1, 5</p> <p>Pulse accuracy ± 5% 1, 5</p> <p><i>Method</i></p> <p>Oscillometric measurement method 1, 5</p> <p>Pulse 40 bpm – 180 bpm 1, 5, 8</p> <p>Manually initiated measurements 13</p> <p>Measurements are from single inflations 13</p> <p><i>Inflation</i></p> <p>Inflation 0 mmHg – 299 mmHg 1, 5, 7</p> <p>Automatic Inflation 7</p> <p>Fuzzy Logic^{Query 2} 7</p> <p>Press button if BP > 220 mmHg 7</p> <p>Manually adjustable inflation pressure 7</p> <p><i>Deflation</i></p> <p>Automatic Deflation 8</p> <p>Automatic safety release valve^{Query 2} 8</p> <p><i>Cuffs</i></p> <p>Single 152 mm × 600 mm (Arm circ. 22 to 42 cm) 6</p> | <p>Measurement</p> <p><i>Accuracy</i></p> <p>BP accuracy ± 3 mmHg^{Query 1} 1, 5</p> <p>Pulse accuracy ± 5% 1, 5</p> <p><i>Method</i></p> <p>Oscillometric measurement method 1, 5</p> <p>Pulse 40 bpm – 180 bpm 1, 5, 8</p> <p>Manually initiated measurements 13</p> <p>Measurements are from single inflations 13</p> <p><i>Inflation</i></p> <p>Inflation 0 mmHg – 299 mmHg 1, 5, 7</p> <p>Automatic Inflation 7</p> <p>Fuzzy Logic^{Query 2} 7</p> <p>Press button if BP > 220 mmHg 7</p> <p>Manually adjustable inflation pressure 7</p> <p><i>Deflation</i></p> <p>Automatic Deflation 8</p> <p>Automatic safety release valve^{Query 2} 8</p> <p><i>Cuffs</i></p> <p>Single 152 mm × 600 mm (Arm circ. 22 to 42 cm) 6</p> |

| Devices | Omron BP785 (HEM-7222-Z) | Omron M6 Comfort (HEM-7000-E) |
|--|--|--|
| Same Criteria (continued) | Buttons/Switches | Buttons/Switches |
| | <i>Power</i> | <i>Power</i> |
| | On/Off with Start/Stop (Start/Stop Label) 10 | On/Off with Start/Stop (O/I Start Label) 10 |
| | <i>Settings</i> | <i>Settings</i> |
| | Date/Time set 10 | Date/Time set 10 |
| | Display/Symbols/Indicators | Display/Symbols/Indicators |
| | <i>Measurement Procedure</i> | <i>Measurement Procedure</i> |
| | Deflation symbol 11 | Deflation symbol 11 |
| | During Measurement: BP Level & Heartbeat 11 | During Measurement: BP Level & Heartbeat 11 |
| | <i>Post Measurement</i> | <i>Post Measurement</i> |
| | SBP, DBP and Pulse 11 | SBP, DBP and Pulse 11 |
| | Average icon 11, 13, 14 | Average icon 11, 13, 14 |
| | Body movement error 3, 11, 13, 18 | Body movement error 3, 11, 13, 18 |
| | Irregular heartbeat 11, 13, 18 | Irregular heartbeat 11, 13, 18 |
| | <i>Measurement Records</i> | <i>Measurement Records</i> |
| | Memory icon 11 | Memory icon 11 |
| | Memory recall number (Replaces pulse rate momentarily) ^{Query 3} 11 | Memory recall number (Replaces pulse rate momentarily) ^{Query 3} 11 |
| | <i>Date and Time</i> | <i>Date and Time</i> |
| | Date and Time 11 | Date and Time 11 |
| | Date and Time (During memory recall) 11 | Date and Time (During memory recall) 11 |
| | <i>Power</i> | <i>Power</i> |
| | Low battery 11, 17 | Low battery 11, 17 |
| | Algorithms | Algorithms |
| | <i>Averages and Differences</i> | <i>Averages and Differences</i> |
| | Last 3 measurements (within 10 min of each other) mean 13 | Last 3 measurements (within 10 min of each other) mean 13 |
| | <i>Diagnostic</i> | <i>Diagnostic</i> |
| | Normotension/Hypertension 13 | Normotension/Hypertension 13 |
| | 135 / 85 mmHg thresholds 13 | 135 / 85 mmHg thresholds 13 |
| Irregular heartbeat detection 13 | Irregular heartbeat detection 13 | |
| Body movement error detection 3, 13 | Body movement error detection 3, 13 | |
| Case | Case | |
| <i>Display</i> | <i>Display</i> | |
| Single screen display 10 | Single screen display 10 | |
| Segment LCD 10 | Segment LCD 10 | |
| <i>Power</i> | <i>Power</i> | |
| 4 "AA" batteries ~ 500 measurements 17 | 4 "AA" batteries ~ 1500 measurements 17 | |
| AC adapter 17 | AC adapter (Optional) 17 | |

| Devices | Omron BP785 (HEM-7222-Z) | Omron M6 Comfort (HEM-7000-E) |
|----------------------------|---|--|
| Comparable Criteria | <p>Measurement Sensors Pressure sensor: piezo-resistive ^{Note 1} 5</p> <p>Measurement Records Memory: 100 measurements × 2 14</p> <p>Buttons/Switches Measurement Records Memory 10</p> <p>Settings Up and down 10</p> <p>Display/Symbols/Indicators Post Measurement Measurement error E_1, E_2, E_3, E_4, E_5 and E_r ^{Query 2, Note 2} 11</p> <p>Hypertension (Indicator strip) 11, 13</p> <p>Case Power Automatic switch-off when not used for 2 min 17</p> | <p>Measurement Sensors Pressure sensor: capacitive ^{Note 1} 5</p> <p>Measurement Records Memory: 90 measurements 14</p> <p>Buttons/Switches Measurement Records Memory × 2 10</p> <p>S U</p> <p>Display/Symbols/Indicators Post Measurement Measurement error $EE/\square, E$ and E/E ^{Query 2, Note 2} 11</p> <p>Hypertension (Blinking heartbeat) 11, 13</p> <p>Case Power Automatic switch-off when not used for 5 min 17</p> |
| Device 2 Criteria | | |


| Queries | | |
|---------|--|--|
| 1 | <p>Query BP accuracy is claimed to be ± 3 mmHg for the M6 Comfort (HEM-7000-E) but ± 3 mmHg or ± 2% for the BP785 (HEM-7222-Z). Can you clarify how these differ, both for pressures above and below 150 mmHg (where the 2% error equals 3 mmHg)?</p> <p>Response <i>BP785 have as same BP accuracy as M6 Comfort (± 3 mmHg). The description “±3mmHg or 2% of reading” comes from the requirement of AAMI SP-10:2008 which is one of the standard for medical device in US. (BP785 are available in US market.) Please find the excerpt from the standard as following.</i></p> <p>4.4.4.B Pressure transducer accuracy At any single condition within the ambient temperature range of 50 °F to 104 °F (10 °C to 40 °C) and the relative humidity range of 15 % to 90 % (non-condensing), both for increasing and for decreasing pressure, the maximum error for the measurement of the cuff pressure at any point of the scale range shall be ± 3 mmHg (± 0.4 kPa) or 2 % of the reading above 200 mmHg.</p> <p>Comment The explanation is accepted</p> | |
| 2 | <p>Query There are differences in the descriptions of the rapid air release, fuzzy logic and error codes between the manuals. Similar queries were raised previously but it is not possible to infer scientifically that answers can be applied in these instances also. Can you confirm that the rapid air release and fuzzy logic are used in both devices and that the mapping of the errors, as described previously for specific devices also applies to these?</p> <p style="text-align: center;">Rapid Air Release Fuzzy Logic Error Codes²</p> | |

| | |
|-------|--|
| | <p>BP785 (HEM-7222-Z) Yes No E1 E4 E5 E2 E3 Er</p> <p>M6 Comfort (HEM-7000-E) No Yes E EE/P E/E</p> <p>Note 1 This is not included in the manual but stated in a previous communication.</p> <p>Note 2 From previous communications, the errors are equivalent and grouped as shown and <i>P</i> refers to a pressure level.</p> <p>Response <i>We confirm that rapid air release function and fuzzy logic are applied for all devices. Regarding to the error codes, please refer the document which we sent previously.</i></p> <p>Comment The explanation is accepted</p> |
| 3 | <p>Query On the respective diagrams for the unit display, the pulse display is also indicated as having a memory number function. This, latter function is not described elsewhere in the manuals. Is this an error in the description or, as described for other devices, does the memory number flash briefly when displaying previous measurements. If so, what is flashed when TruRead averages are shown?</p> <p>Response <i>As described for other devices, the memory number flash briefly when displaying previous measurements. For BP785, memory number is shown before the average is displayed. Memory number is not shown when the values for the individual measurement is displayed.</i></p> <p>Comment The explanation is accepted</p> |
| Notes | <p>1 The sensors used in the Omron M6 Comfort (HEM-7000-E) is known as the “current pressure sensor” (CPSU), a capacitive type, which is being replaced by as new pressure sensor (NPS), a piezoelectric semiconductor type. This is used in the Omron BP785 (HEM-7222-Z).</p> <p>Details of comparatives tests between the sensors have been reviewed by dabl®Educational. Furthermore, the Omron M6 Comfort (HEM-7221-E8), which is the same as the Omron M6 Comfort (HEM-7221-E) except for a similar change in sensor, has been validated using the ESH-IP 2010 protocol and is recommended for use. Following a review of these documents, it was concluded that the change in sensor would not have a detrimental effect on the accuracy of the device.</p> |
| | <p>2 This note from the equivalence application for the HEM-7221-E is also relevant to the HEM-7222-Z.</p> <p><i>Regarding to Group 4, M6 Comfort (7000) error code E had subdivide to M6 Comfort (7221) error code E1, E4 and E5. EE/0 is as same as E2. E/E is as same as E3. The background is explained below. For M6 Comfort (7000), EE/0 is as same as EE, 0 means 0mmHg, and this has the error code Er, but not described in manual. We consider there is no change in the error codes and algorithms among these devices.</i></p> <p><i>For our software, error codes consist of several error judgment conditions. We had a limitation to show enough information on the display in the past due to technical restriction on hardware. For now, the hardware performance has advanced to display more error code. Therefore, we reconsidered the constitution of the error judgment conditions and changed the expression to make it more easy to</i></p> |

understand for users, starting from M6 (HEM-7211-E) and M6 Comfort (HEM-7221-E).

Group 4 Error Codes

| Model | Error codes | | | | | |
|-------------------|-------------|----|-----|----|----|----|
| M6 Comfort (7000) | EE/O | E | E/E | | | |
| M6 Comfort (7221) | E1 | E2 | E3 | E2 | E5 | Er |



| | |
|-----------------------|-----------------------------|
| Recommendation | Equivalence is recommended. |
| Date | 02/07/2012 |